




CCAR1 Antibody

Product Code	CSB-PA816898ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IX12
Immunogen	Recombinant Human Cell division cycle and apoptosis regulator protein 1 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cell division cycle and apoptosis regulator protein 1 (Cell cycle and apoptosis regulatory protein 1) (CARP-1) (Death inducer with SAP domain), CCAR1, CARP1 DIS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CCAR1

Image



Western blot

All lanes: CCAR1 antibody at 8µg/ml

Lane 1: Hela whole cell lysate

Lane 2: PC-3 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 133, 132 kDa

Observed band size: 133 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.