




MICA Antibody

Product Code	CSB-PA013806ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q29983
Immunogen	Recombinant Human MHC class I polypeptide-related sequence A protein (24-307AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	MHC class I polypeptide-related sequence A (MIC-A), MICA, PERB11.1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MICA

Image



Western blot

All lanes: MICA antibody at 8µg/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: A431 whole cell lysate


Lane 3: Mouse Stomach tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 43, 32 kDa

Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded

human colon cancer using CSB-

PA013806ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.