





RSAD2 Antibody

Product Code	CSB-PA848838ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXG1
Immunogen	Recombinant Human Radical S-adenosyl methionine domain-containing protein 2 protein (42-361AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Radical S-adenosyl methionine domain-containing protein 2 (Cytomegalovirus-induced gene 5 protein) (Viperin) (Virus inhibitory protein, endoplasmic reticulum-associated, interferon-inducible), RSAD2, CIG5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RSAD2

Image



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA848838ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA848838ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.