




MSI2 Antibody

Product Code	CSB-PA836197ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DH6
Immunogen	Recombinant Human RNA-binding protein Musashi homolog 2 protein (1-80AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein Musashi homolog 2 (Musashi-2), MSI2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MSI2

Image



Western blot
 All lanes: MSI2 antibody at 10µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 36, 29, 18 kDa
 Observed band size: 36 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.