





SULT1C2 Antibody

Product Code	CSB-PA022938ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00338
Immunogen	Recombinant Human Sulfotransferase 1C2 protein (1-120AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sulfotransferase 1C2 (ST1C2) (EC 2.8.2.-) (Sulfotransferase 1C1) (SULT1C#1) (humSULTC2), SULT1C2, SULT1C1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SULT1C2

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA022938ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA022938ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.