





ATP5O Antibody

Product Code	CSB-PA002379ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48047
Immunogen	Recombinant Human ATP synthase subunit O, mitochondrial protein (24-213AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:1000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP synthase subunit O, mitochondrial (Oligomycin sensitivity conferral protein) (OSCP), ATP5O, ATPO
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ATP5O

Image



Western blot
 All lanes: ATP5O antibody at 3µg/ml + HepG2 whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 23 kDa
 Observed band size: 23 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA002379ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.