





EDA2R Antibody

Product Code	CSB-PA881011ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAV5
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member 27 protein (1-138AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member 27 (X-linked ectodysplasin-A2 receptor) (EDA-A2 receptor), EDA2R, TNFRSF27 XEDAR
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	EDA2R

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA881011ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA881011ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.