





FBXW4 Antibody

Product Code	CSB-PA008525ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57775
Immunogen	Recombinant Human F-box/WD repeat-containing protein 4 protein (143-412AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 4 (Dactylin) (F-box and WD-40 domain-containing protein 4), FBXW4, FBW4 SHFM3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXW4

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA008525ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA008525ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.