




CD6 Antibody

Product Code	CSB-PA326953ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30203
Immunogen	Recombinant Human T-cell differentiation antigen CD6 protein (424-668AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell differentiation antigen CD6 (T12) (TP120) (CD antigen CD6) [Cleaved into: Soluble CD6], CD6
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD6

Image



Immunohistochemistry of paraffin-embedded human thymus tissue using CSB-PA326953ESR2HU at dilution of 1:100