





NLRP2 Antibody

Product Code	CSB-PA865163ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX02
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 2 protein (1-220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 2 (Nucleotide-binding site protein 1) (PYRIN domain and NACHT domain-containing protein 1) (PYRIN-containing APAF1-like protein 2), NLRP2, NALP2 NBS1 PAN1 PYPAF2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NLRP2

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA865163ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA865163ESR2HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.