




TAF5L Antibody

Product Code	CSB-PA023093ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75529
Immunogen	Recombinant Human TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L protein (1-325AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	TAF5-like RNA polymerase II p300/CBP-associated factor-associated factor 65 kDa subunit 5L (PCAF-associated factor 65 beta) (PAF65-beta), TAF5L, PAF65B
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TAF5L

Image



Western blot

All lanes: TAF5L antibody at 10μg/ml + Jurkat

whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 67, 37 kDa

Observed band size: 67 kDa