




L3MBTL1 Antibody

Product Code	CSB-PA897111ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y468
Immunogen	Recombinant Human Lethal(3)malignant brain tumor-like protein 1 protein (513-752AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lethal(3)malignant brain tumor-like protein 1 (H-I(3)mbt) (H-I(3)mbt protein) (L(3)mbt-like) (L(3)mbt protein homolog) (L3MBTL1), L3MBTL1, KIAA0681 L3MBT L3MBTL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	L3MBTL1

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA897111ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.