




Recombinant Mouse Placenta-expressed transcript 1 protein (Plet1)

Product Code	CSB-EP851843MO
Relevance	Modulates leading keratinocyte migration and cellular adhesion to matrix proteins during a wound-healing response and promotes wound repair. May play a role during trichilemmal differentiation of the hair follicle.
Abbreviation	Recombinant Mouse Plet1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VEN2
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SDNGSCVVDNIYTS DILEISTMANVSGGDVTYTVTVPVNDSVSAVILKAVKED DSPVGTSWGTYEKCNDSSVYYNLTSQSQSVFQTNWTVPTSEDVTKVNLQVLI VVNRTASKSSVKMEQVQPSASTPIPESETSQTINTTPTVNTAKTTAKDTANTT AVTTANTTANTTAVTTAKTTAKSLAIRTLGS
Research Area	Cell Biology
Source	E.coli
Target Names	Plet1
Expression Region	28-218aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	36.1kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.