




Recombinant Arabidopsis thaliana Cystathionine beta-lyase, chloroplastic (At3g57050)

Product Code	CSB-EP345319DOA
Relevance	L-cystathionine + H ₂ O = L-homocysteine + NH ₃ + pyruvate.
Abbreviation	Recombinant Mouse-ear cress Cystathionine beta-lyase, chloroplastic protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P53780
Alias	Short name: CBL Alternative name(s): Beta-cystathionase Cysteine lyase
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	NNTTDSLNTMNIKEEASVSTLLVNLDNKFPDAMSTPLYQTATFKQPSAIENG PYDYTRSGNPTRDALESLLAKLDKADRAFCFTSGMAALSAVTHLIKNGEEIVAG DDVYGGSDRLLSQVVPRSGVVVKRVNTTKLDEVAAAIGPQTKLVWLESPTNP RQQISDIRKISEMAHAQGALVLVDNSIMSPVLSRPLELGADIVMHSATKFIAGHS DVMAGVLAVKGEKLAKEVYFLQNSEGSLAPFDCWLCLRGIKTMALRIEKQQE NARKIAMYLSSHPRVKKVYYAGLPDHPGHHLHFSQAKGAGSVFSFITGSVALS KHLVETTKYFSIAVSFGSVKSLISMPCFMASHASIPAEVREARGLTEDLVRISAGI EDVDDLISDLIDIAFKTFPL
Research Area	Others
Source	E.coli
Target Names	At3g57050
Expression Region	70-464aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	46.9kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.