



CA9 Antibody

Product Code	CSB-PA614990GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16790
Immunogen	recombinant Carbonic anhydrase 9 protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,WB,IHC
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	CA9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.