




SLC44A4 Antibody

Product Code	CSB-PA681506LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53GD3
Immunogen	Recombinant Human Choline transporter-like protein 4 protein (56-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Choline transporter-like protein 4 (Solute carrier family 44 member 4) (Thiamine pyrophosphate transporter 1) (hTPPT1), SLC44A4, C6orf29 CTL4 NG22 TPPT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC44A4

Image



Immunofluorescent analysis of HeLa cells using CSB-PA681506LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)