



MACC1 Antibody

Product Code	CSB-PA761627GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZN28
Immunogen	Recombinant Human Metastasis-associated in colon cancer protein 1 protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA,IHC
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	MACC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.