



SMAD4 Antibody

Product Code	CSB-PA619768GA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13485
Immunogen	Recombinant Human Mothers against decapentaplegic homolog 4 protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	SMAD4