



ADAMTS4 Antibody, HRP conjugated

Product Code	CSB-PA001311LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75173
Immunogen	Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 4 protein (213-319AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	A disintegrin and metalloproteinase with thrombospondin motifs 4 (ADAM-TS 4) (ADAM-TS4) (ADAMTS-4) (EC 3.4.24.82) (ADMP-1) (Aggrecanase-1), ADAMTS4, KIAA0688
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ADAMTS4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.