





SERPINB8 Antibody

Product Code	CSB-PA021076EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50452
Immunogen	Recombinant Human Serpin B8 protein (60-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:200-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serpin B8 (Cytoplasmic antiproteinase 2) (CAP-2) (CAP2) (Peptidase inhibitor 8) (PI-8), SERPINB8, PI8
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SERPINB8

Image



Western blot
 All lanes: SERPINB8 antibody at 2µg/ml + Jurkat whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 43, 28, 23 kDa
 Observed band size: 43 kDa



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA021076EA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.