



HSPBAP1 Antibody, FITC conjugated

Product Code	CSB-PA853412LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EW2
Immunogen	Recombinant Human HSPB1-associated protein 1 (1-488AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HSPB1-associated protein 1 (27 kDa heat shock protein-associated protein 1) (Protein associated with small stress protein 1), HSPBAP1, PASS1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HSPBAP1