



HSPBAP1 Antibody, FITC conjugated

Product Code	CSB-PA853412LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EW2
Immunogen	Recombinant Human HSPB1-associated protein 1 (1-488AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	HSPB1-associated protein 1 (27 kDa heat shock protein-associated protein 1) (Protein associated with small stress protein 1), HSPBAP1, PASS1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HSPBAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.