




ACBD4 Antibody

Product Code	CSB-PA818758LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NC06
Immunogen	Recombinant Human Acyl-CoA-binding domain-containing protein 4 protein (1-305AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl-CoA-binding domain-containing protein 4, ACBD4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ACBD4

Image



Western blot

All lanes: ACBD4 antibody at 2µg/ml

Lane 1: Mouse kidney tissue


Lane 2: Mouse liver tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 31, 35, 38 kDa

Observed band size: 31 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA818758LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA818758LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.