



DAO Antibody, HRP conjugated

Product Code	CSB-PA006494LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14920
Immunogen	Recombinant Human D-amino-acid oxidase protein (1-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	D-amino-acid oxidase (DAAO) (DAMOX) (DAO) (EC 1.4.3.3), DAO, DAMOX
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	DAO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.