



ABHD3 Antibody, FITC conjugated

Product Code	CSB-PA845142EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WU67
Immunogen	Recombinant Human Phospholipase ABHD3 protein (52-309AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipase ABHD3 (EC 3.1.1.32) (EC 3.1.1.4) (Abhydrolase domain-containing protein 3), ABHD3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ABHD3