





Cleaved-CASP3 (D175) Antibody

Product Code	CSB-PA000004
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42574
Immunogen	Synthesized peptide derived from the Internal region of Human Caspase-3 p17.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CASP3; CPP32; Caspase-3; CASP-3; Apopain; Cysteine protease CPP32; CPP-32; Protein Yama; SREBP cleavage activity 1; SCA-1
Immunogen Species	Homo sapiens (Human)
Target Names	CASP3

Image



Western Blot analysis of HEK293 cells using Cleaved-Caspase-3 p17 (D175) Polyclonal Antibody



Western Blot analysis of HepG2 cells using Cleaved-Caspase-3 p17 (D175) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.