





Cleaved-CASP6 (D162) Antibody

Product Code	CSB-PA000005
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55212
Immunogen	Synthesized peptide derived from the Internal region of Mouse Caspase-6 p18.
Raised In	Rabbit
Species Reactivity	Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CASP6; MCH2; Caspase-6; CASP-6; Apoptotic protease Mch-2
Immunogen Species	Mus musculus (Mouse)
Target Names	CASP6

Image



Western Blot analysis of HeLa cells using Cleaved-Caspase-6 p18 (D162) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Cleaved-Caspase-6 p18 (D162) Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.