




Cleaved-CASP7 (D198) Antibody

Product Code	CSB-PA000007
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55210
Immunogen	Synthesized peptide derived from the Internal region of Human Caspase-7 p20.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CASP7; MCH3; Caspase-7; CASP-7; Apoptotic protease Mch-3; CMH-1; ICE-like apoptotic protease 3; ICE-LAP3
Immunogen Species	Homo sapiens (Human)
Target Names	CASP7

Image



Western Blot analysis of Jurkat cells using Cleaved-Caspase-7 p20 (D198) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.