




ADAM10 Antibody

Product Code	CSB-PA000013
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14672
Immunogen	Synthesized peptide derived from the Internal region of Human ADAM10.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ADAM10; KUZ; MADM; Disintegrin and metalloproteinase domain-containing protein 10; ADAM 10; CDw156; Kuzbanian protein homolog; Mammalian disintegrin-metalloprotease; CD antigen CD156c
Immunogen Species	Homo sapiens (Human)
Target Names	ADAM10

Image



Western Blot analysis of A549 cells using ADAM10 Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.