





Cleaved-CASP4 (Q81) Antibody

Product Code	CSB-PA000023
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49662
Immunogen	Synthesized peptide derived from the Internal region of Human Caspase-4 p20.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CASP4; ICH2; Caspase-4; CASP-4; ICE(rel)-II; Protease ICH-2; Protease TX
Immunogen Species	Homo sapiens (Human)
Target Names	CASP4

Image



Western Blot analysis of 3T3 cells using Cleaved-Caspase-4 p20 (Q81) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Cleaved-Caspase-4 p20 (Q81) Polyclonal Antibody



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.