




Cleaved-CASP5 (S331) Antibody

Product Code	CSB-PA000024
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51878
Immunogen	Synthesized peptide derived from the C-terminal region of Human Caspase-5 p10.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CASP5; ICH3; Caspase-5; CASP-5; ICE(rel)-III; Protease ICH-3; Protease TY
Immunogen Species	Homo sapiens (Human)
Target Names	CASP5

Image



Western Blot analysis of 293 cells using Cleaved-Caspase-5 p10 (S331) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.