




Cleaved-CTSA (R326) Antibody

Product Code	CSB-PA000028
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10619
Immunogen	Synthesized peptide derived from the Internal region of Human Cathepsin A 32k.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTSA; PPGB; Lysosomal protective protein; Carboxypeptidase C; Carboxypeptidase L; Cathepsin A; Protective protein cathepsin A; PPCA; Protective protein for beta-galactosidase
Immunogen Species	Homo sapiens (Human)
Target Names	CTSA

Image



Western Blot analysis of Jurkat cells using Cleaved-Cathepsin A 32k (R326) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.