




Cleaved-CTSC (R394) Antibody

Product Code	CSB-PA000029
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53634
Immunogen	Synthesized peptide derived from the C-terminal region of Human Cathepsin C HC.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTSC; CPPI; Dipeptidyl peptidase 1; Cathepsin C; Cathepsin J; Dipeptidyl peptidase I; DPP-I; DPPI; Dipeptidyl transferase
Immunogen Species	Homo sapiens (Human)
Target Names	CTSC

Image



Western Blot analysis of Jurkat cells using Cleaved-Cathepsin C HC (R394) Polyclonal Antibody