




Cleaved-CTSL (T288) Antibody

Product Code	CSB-PA000034
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07711
Immunogen	Synthesized peptide derived from the C-terminal region of Human Cathepsin L1 HC.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTSL1; CTSL; Cathepsin L1; Major excreted protein; MEP
Immunogen Species	Homo sapiens (Human)
Target Names	CTSL

Image



Western Blot analysis of HeLa cells using Cleaved-Cathepsin L1 HC (T288) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.