




Cleaved-CTSZ (L62) Antibody

Product Code	CSB-PA000038
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBR2
Immunogen	Synthesized peptide derived from the Internal region of Human Cathepsin Z.
Raised In	Rabbit
Species Reactivity	Human,Monkey
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CTSZ; Cathepsin Z; Cathepsin P; Cathepsin X
Immunogen Species	Homo sapiens (Human)
Target Names	CTSZ

Image



Western Blot analysis of COS7 cells using Cleaved-Cathepsin Z (L62) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.