




DFFA Antibody

Product Code	CSB-PA000044
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00273
Immunogen	Synthesized peptide derived from the Internal region of Human ICAD.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	DFFA; DFF1; DFF45; H13; DNA fragmentation factor subunit alpha; DNA fragmentation factor 45 kDa subunit; DFF-45; Inhibitor of CAD; ICAD
Immunogen Species	Homo sapiens (Human)
Target Names	DFFA

Image



Western Blot analysis of 293 cells using ICAD Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.