





Cleaved-MMP10 (F99) Antibody

Product Code	CSB-PA000048
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09238
Immunogen	Synthesized peptide derived from the Internal region of Human MMP-10.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:40000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MMP10; STMY2; Stromelysin-2; SL-2; Matrix metalloproteinase-10; MMP-10; Transin-2
Immunogen Species	Homo sapiens (Human)
Target Names	MMP10

Image



Western Blot analysis of HepG2 cells using Cleaved-MMP-10 (F99) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Cleaved-MMP-10 (F99) Polyclonal Antibody