




Cleaved-MMP23A (Y79) Antibody

Product Code	CSB-PA000054
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75900
Immunogen	Synthesized peptide derived from the Internal region of Human MMP-23.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	MMP23A; MMP21; MMP23B; MMP21; MMP22; Matrix metalloproteinase-23; MMP-23; Femalysin; MIFR-1; Matrix metalloproteinase-21; MMP-21; Matrix metalloproteinase-22; MMP-22
Immunogen Species	Homo sapiens (Human)
Target Names	MMP23A

Image



Western Blot analysis of 293 cells using Cleaved-MMP-23 (Y79) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.