





Cleaved-NOTCH1 (V1754) Antibody

Product Code	CSB-PA000057
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46531
Immunogen	Synthesized peptide derived from the Internal region of Human Notch 1.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	NOTCH1; TAN1; Neurogenic locus notch homolog protein 1; Notch 1; hN1; Translocation-associated notch protein TAN-1
Immunogen Species	Homo sapiens (Human)
Target Names	NOTCH1

Image



Western Blot analysis of KB cells using Cleaved-Notch 1 (V1754) Polyclonal Antibody



Western Blot analysis of MOUSE-BRAIN RAT-MUSCLE cells using Cleaved-Notch 1 (V1754) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.