




Cleaved-F7 (R212) Antibody

Product Code	CSB-PA000066
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08709
Immunogen	Synthesized peptide derived from the Internal region of Human Factor VII LC.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	F7; Coagulation factor VII; Proconvertin; Serum prothrombin conversion accelerator; SPCA; Eptacog alfa
Immunogen Species	Homo sapiens (Human)
Target Names	F7

Image



Western Blot analysis of Jurkat cells using Cleaved-Factor VII LC (R212) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.