




HABP2 Antibody

Product Code	CSB-PA000069
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14520
Immunogen	Synthesized peptide derived from the Internal region of Human HABP2.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	HABP2; HGFAL; PHBP; Hyaluronan-binding protein 2; Factor VII-activating protease; Factor seven-activating protease; FSAP; Hepatocyte growth factor activator-like protein; Plasma hyaluronan-binding protein
Immunogen Species	Homo sapiens (Human)
Target Names	HABP2

Image



Western Blot analysis of HeLa cells using HABP2 Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.