




# Cleaved-ITGA6 (E942) Antibody

<b>Product Code</b>	CSB-PA000072
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P23229
<b>Immunogen</b>	Synthesized peptide derived from the C-terminal region of Human Integrin $\alpha 6$ LC.
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Purification Method</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Isotype</b>	IgG
<b>Alias</b>	ITGA6; Integrin alpha-6; CD49 antigen-like family member F; VLA-6; CD antigen CD49f
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ITGA6

## Image



Western Blot analysis of HeLa cells using Cleaved-Integrin  $\alpha 6$  LC (E942) Polyclonal Antibody

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.