




Cleaved-ITGAV (K889) Antibody

Product Code	CSB-PA000074
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06756
Immunogen	Synthesized peptide derived from the C-terminal region of Human Integrin α V HC.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ITGAV; MSK8; VNRA; Integrin alpha-V; Vitronectin receptor subunit alpha; CD antigen CD51
Immunogen Species	Homo sapiens (Human)
Target Names	ITGAV

Image



Western Blot analysis of A549 cells using Cleaved-Integrin α V HC (K889) Polyclonal Antibody