





HA-tag Antibody

Product Code	CSB-PA000087
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	HA synthetic peptide conjugated to KLH.
Raised In	Rabbit
Tested Applications	WB, IP, ELISA; WB:1:1000-1:3000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Product Type	Polyclonal Antibody

Image



Western Blot analysis of HA-protein cells using HA-tag Polyclonal Antibody



IP analysis of HA-tag Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.