





Acetyl-Histone H2B (K5) Antibody

Product Code	CSB-PA000098
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57053
Immunogen	Synthesized peptide derived from Human Histone H2B around the acetylation site of K5.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	H2BFS; Histone H2B type F-S; Histone H2B.s; H2B/s
Immunogen Species	Homo sapiens (Human)
Target Names	H2BFS

Image



Western Blot analysis of 3T3 cells using Acetyl-Histone H2B (K5) Polyclonal Antibody



Western Blot analysis of VEC cells using Acetyl-Histone H2B (K5) Polyclonal Antibody