





Acetyl-Histone H3 (K23) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA000100 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P68431/Q71DI3/P84243 |
| Immunogen | Synthesized peptide derived from Human Histone H3 around the acetylation site of K23. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, IF, ELISA; WB:1:500-1:2000, IF:1:200-1:1000, ELISA:1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | HIST1H3A; H3FA; HIST1H3B; H3FL; HIST1H3C; H3FC; HIST1H3D; H3FB; HIST1H3E; H3FD; HIST1H3F; H3FI; HIST1H3G; H3FH; HIST1H3H; H3FK; HIST1H3I; H3FF; HIST1H3J; H3FJ; Histone H3.1; Histone H3/a; Histone H3/b; Histone H3/c; Histone H3/d; Histone H3 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HIST1H3A/HIST2H3A/H3F3A |

Image



Western Blot analysis of L929 cells using Acetyl-Histone H3 (K23) Polyclonal Antibody



Western Blot analysis of L929 cells using Acetyl-Histone H3 (K23) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.