





Acetyl-TP73 (K327) Antibody

Product Code	CSB-PA000126
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15350
Immunogen	Synthesized peptide derived from the Human p73 around the acetylation site of K327.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP73; P73; Tumor protein p73; p53-like transcription factor; p53-related protein
Immunogen Species	Homo sapiens (Human)
Target Names	TP73

Image



Western Blot analysis of JK cells using Acetyl-p73 (K327) Polyclonal Antibody



Western Blot analysis of Jurkat cells using
Acetyl-p73 (K327) Polyclonal Antibody