




Acetyl-ARNTL (K538) Antibody

Product Code	CSB-PA000129
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00327
Immunogen	Synthesized peptide derived from Human BMAL1 around the acetylation site of K538.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	ARNTL; BHLHE5; BMAL1; MOP3; PASD3; Aryl hydrocarbon receptor nuclear translocator-like protein 1; Basic-helix-loop-helix-PAS protein MOP3
Immunogen Species	Homo sapiens (Human)
Target Names	ARNTL

Image



Western Blot analysis of 293 cells using Acetyl-BMAL1 (K538) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.