





Acetyl-TP53 (K382) Antibody

Product Code	CSB-PA000136
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Synthesized peptide derived from the Human p53 around the acetylation site of K382.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	TP53; P53; Cellular tumor antigen p53; Antigen NY-CO-13; Phosphoprotein p53; Tumor suppressor p53
Immunogen Species	Homo sapiens (Human)
Target Names	TP53

Image



Western Blot analysis of 293T-UV cells using Acetyl-p53 (K382) Polyclonal Antibody



Western Blot analysis of HepG2 cells using Acetyl-p53 (K382) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.