





Acetyl-UBA52/RPS27A/UBB/UBC (K27) Antibody

Product Code	CSB-PA000137
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62987/P62979/P0CG47/P0CG48
Immunogen	Synthesized peptide derived from the Human Ub around the acetylation site of K27.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	UBB; Polyubiquitin-B; UBC; Polyubiquitin-C; RPS27A; UBA80; UBCEP1; Ubiquitin-40S ribosomal protein S27a; Ubiquitin carboxyl extension protein 80; UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion pr
Immunogen Species	Homo sapiens (Human)
Target Names	UBA52/RPS27A/UBB/UBC

Image



Western Blot analysis of HEPG2 cells using Acetyl-Ub (K27) Polyclonal Antibody



Western Blot analysis of HeLa cells using Acetyl-Ub (K27) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.